


Cranial Plate Set II

iSpine

The iSpine® Cranial plating system is a low-profile line of plates, screws and meshes intended for craniotomy and rigid fixation of craniofacial fractures.



Low Profile Fixation System

- *High torque and low-profile fixation screws*
- *Wide selection of fixation plates & mesh plates*
- *Ergonomic designed complete set of instruments*

Aesthetics


- *Compact & simple designs of all screws and plates*
- *Smooth and rounded plate edges*

Reliability

- *Slim and strong products*
- *Stable and effective system*
- *Standardised instruments*
- *Reduced inventory for hospitals*

Comfort

- *Easy, safe and accurate insertion*
- *Precision screws clip onto screw driver head and do not drop*
- *Decreased profile, rounded edges & low-profile head screws minimises plate palpability whilst providing strength for cranial fixation*



Bone Plate

- Various types and strengths for optimum operation
- Smooth curves with rounded edges
- High quality fixation plates
- Pure titanium - minimal MRI artifact

Bone Screw

- Self-drilling
- No screw dropping
- Precision fixation
- Stable and easy insertion
- Cross drive
- Emergency screw available
- Titanium alloy

Disposable Battery Powered Driver

- Battery-powered screwdriver in sterile package
- Single-use disposable design eliminates need for autoclaving and repairing
- Designed to drive self-drilling screw

Compatibility

Driver Shaft	Pilot Drill
Micro (1.4mm)	Micro (1.0mm)
Mid (1.5mm)	Mid (1.3mm)


Specifications

Torque at the motor: 30 N·cm
Motor speed: 160 rpm

Each set is equipped with:


1x BD-05 small grasping forceps	4x SDB-15c Screwdriver Shaft Long	32x 4mm screws
1x DP-96 Pliers	2x ADB-15c Screwdriver Shaft Short	36x 5mm screws
1x DP-91 scissors	10 units of each plate	20x 6mm screws
1x 112 MD STOP drill bit	1x 3D Mesh	4x 4mm Rescue screws
1x ADB disposable battery operated driver	1x screen Mesh	4x 5mm Rescue Screws
2x SDH4 screwdriver handles	20x 3mm screws	4x 6mm Rescue Screws

Instrument Log




Self-Drilling Fixation Screw - 1.5 x 3mm
 Self-Drilling Fixation Screw - 1.5 x 4mm
 Self-Drilling Fixation Screw - 1.5 x 5mm
 Self-Drilling Fixation Screw - 1.5 x 6mm

Self-Drilling Fixation Screw - Emergency 1.8 x 4mm
 Self-Drilling Fixation Screw - Emergency 1.8 x 5mm
 Self-Drilling Fixation Screw - Emergency 1.8 x 6mm



1

1. Screen Mesh



2

2. 3D Mesh




1. Plate - 1.5mm x 2 hole straight mid
 2. Plate - 1.5mm x 2 hole straight large
 3. Plate - 1.5mm x 4 hole straight

4. Plate - 1.5mm x 6 hole straight
 5. Plate - 1.5mm x 6 hole double-Y




1




2

1. Plate - 1.5mm x 4 hole X Type
 2. Plate - 1.5mm x 4 hole Square



3


3. Plate - 1.5mm x 6 hole double square
 4. Plate - 1.5mm x 8 hole triple square




4



1



2




3


1. Plate - 1.5mm x 5 hole Burr 12mm

2. Plate - 1.5mm x 6 hole Burr 12mm

3. Plate - 1.5mm x 5 hole Burr 17mm



1



2



3

1. Plate - 1.5mm x 6 hole Burr 17mm

2. Plate - 1.5mm x 6 hole Burr 22mm

3. Plate - 1.5mm x 5 hole Burr 22mm



Disclaimer

This document is intended exclusively for physicians and is not intended for laypersons.

Information on the products and procedures contained in this document is of a general nature and does not represent nor constitute medical advice or recommendations. Because this information does not purport to constitute a diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and / or advice in whole or in part.

Information in this document was gathered and compiled by medical experts and qualified iSpine personnel. The information contained herein is accurate to the best knowledge of iSpine and of those experts and personnel involved in its compilation. However, iSpine does not assume any liability for the accuracy, completeness or quality of the information in this document, and iSpine is not liable for any losses, tangible or intangible, that may be caused by the use of this information.

www.ispine.co.za • admin@ispine.co.za • Mulberry Hill Office Park, Unit 10, Broadacres Drive, Dainfern, Midrand
+27 11 469 0060

24 Hour Emergency Number: +27 82 881 8929